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Orbitrap Exploris 120 质谱仪能够在控制软件 (Xcalibur, 版本：4.4, Thermo) 控制下进行一级、二级质谱数据采集。详细参数如下：sheath gas flow rate as 50 Arb, Aux gas flow rate as 15 Arb, capillary temperature 320 °C, full MS resolution as 60000, MS/MS resolution as 15000, collision energy: SNCE 20/30/40, spray voltage as 3.8 kV (positive) or -3.4 kV (negative), respectively.

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#### 1.Metabolites Extraction:

##### 1.1.Metabolites Extraction of Liquid:

100  $\mu$ L of sample was taken, mixed with 400  $\mu$ L of extraction solution (MeOH:ACN, 1:1 (v/v)), the extraction solution contain deuterated internal standards, the mixed solution were vortexed for 30 s, sonicated for 10 min in 4  $^{\circ}$ C water bath, and incubated for 1 h at -40  $^{\circ}$ C to precipitate proteins. Then the samples were centrifuged at 12000 rpm (RCF=13800( $\times$ g), R= 8.6cm) for 15 min at 4  $^{\circ}$ C. The supernatant was transferred to a fresh glass vial for analysis. The quality control (QC) sample was prepared by mixing an equal aliquot of the supernatant of samples.

##### 1.2 Metabolites Extraction of Solid:

The animal tissue samples (25 mg $\pm$ 1 mg) were taken, mixed with beads and 500  $\mu$ L of extraction

solution (MeOH:ACN:H<sub>2</sub>O, 2:2:1 (v/v)). The extraction solution contained deuterated internal standards. The mixed solution was vortexed for 30 s.

The soil samples (100 mg ± 1 mg) were taken, mixed with beads and 500 µL of extraction solution (MeOH:ACN:H<sub>2</sub>O, 2:2:1 (v/v)). The extraction solution contained deuterated internal standards. The mixed solution was vortexed for 30 s.

The plant samples (20 mg ± 1 mg) were taken and lyophilized, mixed with beads and 1000 µL of extraction solution (MeOH:ACN:H<sub>2</sub>O, 2:2:1 (v/v)). The extraction solution contained deuterated internal standards. The mixed solution was vortexed for 30 s.

Then the mixed samples were homogenized (35 Hz, 4 min) and sonicated for 5 min in 4 °C water bath, the step repeated for three times.

The samples were incubated for 1 h at -40 °C to precipitate proteins. Then the samples were centrifuged at 12000 rpm (RCF=13800(×g), R= 8.6cm) for 15 min at 4 °C. The supernatant was transferred to a fresh glass vial for analysis. The quality control (QC) sample was prepared by mixing an equal aliquot of the supernatant of samples.

### 1.3 Metabolites Extraction of Cell or Bacteria:

The cell pellets (about 10<sup>7</sup> cells) or bacteria pellets (about 10<sup>7</sup> bacteria) were taken, mixed with 1000 µL of extraction solution (MeOH:ACN:H<sub>2</sub>O, 2:2:1 (v/v)), the extraction solution contained deuterated internal standards, the mixed solution was vortexed for 30 s and incubated in liquid nitrogen for 1 min. The samples were then allowed to thaw at room temperature and vortexed for 30 s. This freeze - thaw cycle was repeated three times. Then the samples were sonicated for 10 min in 4 °C water bath, and incubated for 1 h at -40 °C to precipitate proteins. The samples were centrifuged at 12000 rpm (RCF=13800(×g), R= 8.6cm) for 15 min at 4 °C. The supernatant was transferred to a fresh glass vial for analysis. The quality control (QC) sample was prepared by mixing an equal aliquot of the supernatant of samples.

## 2. LC-MS/MS Analysis:

LC-MS/MS analyses were performed using an UHPLC system (Vanquish, Thermo FisherScientific) with a Waters ACQUITY UPLC BEH Amide (2.1 mm × 50 mm, 1.7 μm) coupled to Orbitrap Exploris 120 mass spectrometer (Orbitrap MS, Thermo). The mobile phase consisted of 25 mmol/L ammonium acetate and 25 mmol/L ammonia hydroxide in water (pH = 9.75) (A) and acetonitrile (B). The auto-sampler temperature was 4 °C, and the injection volume was 2 μL. The Orbitrap Exploris 120 mass spectrometer was used for its ability to acquire MS/MS spectra on information-dependent acquisition (IDA) mode in the control of the acquisition software (Xcalibur, Thermo). In this mode, the acquisition software continuously evaluates the full scan MS spectrum. The ESI source conditions were set as following: sheath gas flow rate as 50 Arb, Aux gas flow rate as 15 Arb, capillary temperature 320 °C, full MS resolution as 60000, MS/MS resolution as 15000, collision energy: SNCE 20/30/40, spray voltage as 3.8 kV (positive) or -3.4 kV (negative), respectively.

### 4.1.2. Data preprocessing and annotation:

The raw data were converted to the mzXML format using ProteoWizard and processed with an in-house program, which was developed using R and based on The raw data were converted to the mzXML format using ProteoWizard and processed with an in-house program, which was developed using R and based on XCMS, for peak detection, extraction, alignment, and integration. Then an in-house MS2 database was applied in metabolite annotation. The cutoff for annotation was set at 0.3.

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The acquired MS data pretreatments including peak picking, peak grouping, retention time correction, second peak grouping, and annotation of isotopes and adducts was performed using XCMS software. LC-MS raw data files were converted into mzXML format and then processed by the XCMS, CAMERA and metaX toolbox implemented with the R software. Each ion was identified by combining retention time (RT) and m/z data. Intensities of each peaks were recorded and a three dimensional matrix containing arbitrarily assigned peak indices (retention time-m/z pairs), sample names (observations) and ion intensity information (variables) was generated.

The online KEGG, HMDB database was used to annotate the metabolites by matching the exact molecular mass data (m/z) of samples with those from database. If a mass difference between observed and the database value was less than 10 ppm, the metabolite would be annotated and the molecular formula of metabolites would further be identified and validated by the isotopic distribution measurements. We also used a in-house fragment spectrum library of metabolites to validate the metabolite identification.

Statistical analysis was performed in R (version 4.0.0). The raw protein intensity will be normalized by method "medium", Hierarchical clustering was performed using pheatmap package. Principal component analysis (PCA) was performed using metaX package. The PLSDA analysis is performed by the R package ropls and the VIP values of each variable are calculated. Correlation analysis was performed by Pearson correlation coefficient of cor package. The three conditions of P Value < 0.05, difference multiple > 1.2 obtained by T test and VIP calculated by PLSDA analysis simultaneously met the screening of the final metabolites with significant differences..

Hypergeometric-based enrichment analysis with KEGG Pathway was performed to annotate protein sequences. individually. The software GSEA (v4.1.0) and MSigDB were used for gene set enrichment analysis to determine whether a set of genes in a specific KEGG pathway in different situations. Meeting this condition  $|NES| > 1$ , NOM p-val < 0.05, FDR q-val < 0.25 were considered to be significantly different between the two groups. The network map is drawn according to the pathway where the metabolite is located.

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